


Co-learning and Self-supervised Learning

[Spring 2020 CS-8395 Deep Learning in Medical Image Computing]


Instructor: Yuankai Huo, Ph.D.
Department of Electrical Engineering and Computer Science
Vanderbilt University

Simple Recipe for Classification

“CONV-POOL-FC”




Fully Convolutional Network (FCN)



Encoder


Decoder

FCN:





Long, et al. CVPR (2015)

Unet:



Ronneberger, et al. MICCAI (2015)

Feed-forward and Recurrent Network



<https://towardsdatascience.com/recurrent-neural-networks-and-lstm-4b601dd822a5>

<https://www.leiphone.com/news/201608/syAwLNx4bGPuFYI1.html> <https://slideplayer.com/slide/3383596/>

RNN


$$\mathbf{h}_t = \theta \phi(\mathbf{h}_{t-1}) + \theta_x \mathbf{x}_t$$

$$\mathbf{y}_t = \theta_y \phi(\mathbf{h}_t)$$


$$E = \dots + E_{t-1} + E_t + E_{t+1} + E_{t+2} + \dots$$

$$\frac{\partial E}{\partial \theta} = \sum_{t=1}^S \frac{\partial E_t}{\partial \theta}$$



$$\frac{\partial E_t}{\partial \theta} = \sum_{k=1}^t \frac{\partial E_t}{\partial \mathbf{y}_t} \frac{\partial \mathbf{y}_t}{\partial \mathbf{h}_t} \frac{\partial \mathbf{h}_t}{\partial \mathbf{h}_k} \frac{\partial \mathbf{h}_k}{\partial \theta}$$





RNN




LSTM



LSTM



LSTM



$$f_t = \sigma(W_f \cdot [h_{t-1}, x_t] + b_f)$$


$$i_t = \sigma(W_i \cdot [h_{t-1}, x_t] + b_i)$$

$$\tilde{C}_t = \tanh(W_C \cdot [h_{t-1}, x_t] + b_C)$$

$$C_t = f_t * C_{t-1} + i_t * \tilde{C}_t$$



$$o_t = \sigma(W_o \cdot [h_{t-1}, x_t] + b_o)$$

$$h_t = o_t * \tanh(C_t)$$





<http://colah.github.io/posts/2015-08-Understanding-LSTMs/>


Attention Model




Attention Model with LSTM




Class Activation Mapping




Knowledge For Today's Class




CNN



RNN (LSTM)




CNN with Attention



RNN with Attention

CS-8395 Deep Learning in Medical Image Computing

What do doctor do?



<https://gfycat.com/uk/gifs/search/dr+z>

In Practice



<https://thealthtech.com/oscar-emr/>

How?



<https://www.cdc.gov/features/medical-imaging-procedures/index.html>

<https://www.healthcare-informatics.com/news-item/report-emr-market-hit-25b-2014>

Summarize



CONFIDENTIAL																							
VA Medication History																							
Source: VA																							
Last Updated: 18 Jul 2012 @ 00:00																							
Sorted by: Last Filled On (Descending)																							
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.																							
<table border="1"><thead><tr><th>Medication:</th><th>DOCGATAN 100MG CAP</th></tr></thead><tbody><tr><td>Instructions:</td><td>TAKE ONE CAPSULE BY MOUTH EVERY DAY</td></tr><tr><td>Status:</td><td>Active</td></tr><tr><td>Refills Remaining:</td><td>3</td></tr><tr><td>Last Filled On:</td><td>12 Apr 2012</td></tr><tr><td>Initially Ordered On:</td><td>12 Apr 2012</td></tr><tr><td>Quantity</td><td>Days Supply</td><td>Pharmacy</td><td>Prescription Number</td></tr><tr><td>30</td><td>30</td><td>PORTLAND PHARMACY</td><td>112024128</td></tr></tbody></table>				Medication:	DOCGATAN 100MG CAP	Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Status:	Active	Refills Remaining:	3	Last Filled On:	12 Apr 2012	Initially Ordered On:	12 Apr 2012	Quantity	Days Supply	Pharmacy	Prescription Number	30	30	PORTLAND PHARMACY	112024128
Medication:	DOCGATAN 100MG CAP																						
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY																						
Status:	Active																						
Refills Remaining:	3																						
Last Filled On:	12 Apr 2012																						
Initially Ordered On:	12 Apr 2012																						
Quantity	Days Supply	Pharmacy	Prescription Number																				
30	30	PORTLAND PHARMACY	112024128																				
<table border="1"><thead><tr><th>Medication:</th><th>METOPROLOL TARTRATE 50MG TAB</th></tr></thead><tbody><tr><td>Instructions:</td><td>TAKE ONE TABLET BY MOUTH EVER 12 HOURS</td></tr><tr><td>Status:</td><td>Active</td></tr><tr><td>Refills Remaining:</td><td>3</td></tr><tr><td>Last Filled On:</td><td>17 Apr 2012</td></tr></tbody></table>				Medication:	METOPROLOL TARTRATE 50MG TAB	Instructions:	TAKE ONE TABLET BY MOUTH EVER 12 HOURS	Status:	Active	Refills Remaining:	3	Last Filled On:	17 Apr 2012										
Medication:	METOPROLOL TARTRATE 50MG TAB																						
Instructions:	TAKE ONE TABLET BY MOUTH EVER 12 HOURS																						
Status:	Active																						
Refills Remaining:	3																						
Last Filled On:	17 Apr 2012																						



Summarize



VA Medication History			
Source: VA			
Last Updated: 18 Jul 2012 @ 00:00			
Sort By: Last Filled On (Descending)			
VA Medication history includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCGATE NA 100MG CAP	Days Supply:	30
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11303413A
Refills Remaining:	10	Last Filled On:	12 Apr 2012
Initially Ordered On:	12 Apr 2012	Quantity:	30
Medication:	METOPROLOL TARTRATE 50MG TAB	Days Supply:	30
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11303413A
Refills Remaining:	3	Last Filled On:	12 Apr 2012

VA Medication History			
Source: VA			
Last Updated: 18 Jul 2012 @ 00:00			
Sort By: Last Filled On (Descending)			
VA Medication history includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCGATE NA 100MG CAP	Days Supply:	30
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11303413A
Refills Remaining:	10	Last Filled On:	12 Apr 2012
Initially Ordered On:	12 Apr 2012	Quantity:	30
Medication:	METOPROLOL TARTRATE 50MG TAB	Days Supply:	30
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11303413A
Refills Remaining:	3	Last Filled On:	12 Apr 2012

Summarize




VA Medication History								
Source: VA								
Last Updated: 18 Jul 2012 @ 0000								
Sorted By: Last Filled On (Descending)								
VA Medication history includes up to two years of medication history unless you select a different date range in your download request.								
Medication: DOCGATE NA 100MG CAP								
Instructions: TAKE ONE CAPSULE BY MOUTH EVERY DAY								
Status: Active								
Refills Remaining: 10								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					
Medication: METOPROLOL TARTRATE 50MG TAB								
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS								
Status: Active								
Refills Remaining: 3								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					

VA Medication History								
Source: VA								
Last Updated: 18 Jul 2012 @ 0000								
Sorted By: Last Filled On (Descending)								
VA Medication history includes up to two years of medication history unless you select a different date range in your download request.								
Medication: DOCGATE NA 100MG CAP								
Instructions: TAKE ONE CAPSULE BY MOUTH EVERY DAY								
Status: Active								
Refills Remaining: 10								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					
Medication: METOPROLOL TARTRATE 50MG TAB								
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS								
Status: Active								
Refills Remaining: 3								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					



VA Medication History								
Source: VA								
Last Updated: 18 Jul 2012 @ 0000								
Sorted By: Last Filled On (Descending)								
VA Medication history includes up to two years of medication history unless you select a different date range in your download request.								
Medication: DOCGATE NA 100MG CAP								
Instructions: TAKE ONE CAPSULE BY MOUTH EVERY DAY								
Status: Active								
Refills Remaining: 10								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					
Medication: METOPROLOL TARTRATE 50MG TAB								
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS								
Status: Active								
Refills Remaining: 3								
Last Filled On: 12 Apr 2012								
Initially Ordered On: 12 Apr 2012								
Quantity	Days Supply	Pharmacy	Prescription Number					
30	30	PORTLAND PHARMACY	1130343A					

Strategy



MHVZZOED, TEST	CONFIDENTIAL	Page 11 of 27																				
VA Medication History																						
Source: VA																						
Last Updated: 18 Jul 2012 @ 0002																						
Sorted By: Last Filled On [Descending]																						
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.																						
<table border="1"><tr><td>Medication:</td><td>DOCUSTATE NA 100MG CAP</td></tr><tr><td>Instructions:</td><td>TAKE ONE CAPSULE BY MOUTH EVERY DAY</td></tr><tr><td>Status:</td><td>Active</td></tr><tr><td>Refills Remaining:</td><td>10</td></tr><tr><td>Last Filled On:</td><td>12 Apr 2012</td></tr><tr><td>Initially Ordered On:</td><td>12 Apr 2012</td></tr><tr><td>Quantity</td><td>Days Supply</td><td>Pharmacy</td><td>Prescription Number</td></tr><tr><td>30</td><td>30</td><td>PORTLAND PHARMACY</td><td>11305543A</td></tr></table>			Medication:	DOCUSTATE NA 100MG CAP	Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Status:	Active	Refills Remaining:	10	Last Filled On:	12 Apr 2012	Initially Ordered On:	12 Apr 2012	Quantity	Days Supply	Pharmacy	Prescription Number	30	30	PORTLAND PHARMACY	11305543A
Medication:	DOCUSTATE NA 100MG CAP																					
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY																					
Status:	Active																					
Refills Remaining:	10																					
Last Filled On:	12 Apr 2012																					
Initially Ordered On:	12 Apr 2012																					
Quantity	Days Supply	Pharmacy	Prescription Number																			
30	30	PORTLAND PHARMACY	11305543A																			
<table border="1"><tr><td>Medication:</td><td>METOPROLOL TARTRATE 50MG TAB</td></tr><tr><td>Instructions:</td><td>TAKE ONE TABLET BY MOUTH EVERY 12 HOURS</td></tr><tr><td>Status:</td><td>Active</td></tr><tr><td>Refills Remaining:</td><td>3</td></tr></table>			Medication:	METOPROLOL TARTRATE 50MG TAB	Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS	Status:	Active	Refills Remaining:	3												
Medication:	METOPROLOL TARTRATE 50MG TAB																					
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS																					
Status:	Active																					
Refills Remaining:	3																					
Last Filled On: 12 Apr 2012																						

Strategy



CONFIDENTIAL


Page 11 of 27

VA Medication History



Source: VA
Last Updated: 18 Jul 2012 @ 0002
Sorted By: Last Filled On [Descending]
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.

Medication:	DOCUCLATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12 Apr 2012		
Initially Ordered On:	12 Apr 2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A

Medication:	METOPROLOL TARTRATE 50MG TAB
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS
Status:	Active
Refills Remaining:	3



Strategy






CONFIDENTIAL

Page 11 of 27

VA Medication History			
Source: VA			
Last Updated: 18 Jul 2012 @ 0002			
Sorted By: Last Filled On [Descending]			
Medication:	DOSUCATE NA 100MG CAP	Days Supply:	30
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11305543A
Refills Remaining:	10	Last Filled On:	12 Apr 2012
Initially Ordered On:	12 Apr 2012		
Quantity	30	Days Supply	30
Medication: METOPROLOL TARTRATE 50MG TAB			
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS			
Status: Active			
Refills Remaining: 3			


Strategy



Cancer?

VA Medication History			
Source:	VA		
Last Updated:	18-Jul 2012 @ 0002		
Sorted By:	Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOSUCATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12-Apr-2012		
Initially Ordered On:	12-Apr-2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A
Medication:	METOPROLOL TARTRATE 50MG TAB		
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status:	Active		
Refills Remaining:	3		

Strategy



↓ CNN



↓ DNN or RNN

Cancer?

VA Medication History			
Source:	VA		
Last Updated:	18-Jul 2012 @ 0002		
Sorted By:	Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCUSATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12-Apr-2012		
Initially Ordered On:	12-Apr-2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A
Medication:	METOPROLOL TARTRATE 50MG TAB		
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status:	Active		
Refills Remaining:	3		

CONFIDENTIAL

Page 11 of 27

↓ DNN or RNN



↓ DNN or RNN

Cancer?

Image + Text

TandemNet: Distilling Knowledge from Medical
Images Using Diagnostic Reports
as Optional Semantic References

Zizhao Zhang, Pingjun Chen, Manish Sapkota, Lin Yang

University of Florida




Fig. 1: The illustration of the TandemNet.

Image + Text

Method	Accuracy (%)	
	w/o text	w/ text
WRN16-4	75.4	-
ResNet18-TL	79.4	-
TandemNet-WVS	79.4	85.6
TandemNet	82.4	89.9
TandemNet-TL	84.9	88.6

Table 1: The quantitative evaluation (averaged on 3 trials). The first block shows standard CNNs so text is irrelevant.





Fig. 2: The confusion matrices of two compared methods ResNet18-TL and TandemNet-TL (w/o text) in Table 1.

Strategy



↓ CNN



With Attention



↓ DNN or RNN
Cancer?

MHVZ0ED, TEST	CONFIDENTIAL	Page 11 of 27
<i>VA Medication History</i>		
Source: VA		
Last Updated: 18 Jul 2012 @ 0002		
Sorted By: Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.		
Medication: DOCUSATE NA 100MG CAP		
Instructions: TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status: Active		
Refills Remaining: 10		
Last Filled On: 12 Apr 2012		
Initially Ordered On: 12 Apr 2012		
Quantity	Days Supply	Pharmacy
30	30	PORTLAND PHARMACY
Prescription Number: 11305543A		
Medication: METOPROLOL TARTRATE 50MG TAB		
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status: Active		
Refills Remaining: 3		
Last Filled On: 12 Apr 2012		

↓ DNN or RNN



Be Careful

Are we fooling ourselves?

- E.g. action recognition
 - Very hard to improve on single frame classifiers
 - Consider "opening fridge" action.

example by David Fouhey

A video frame showing a man in a blue sweater standing and speaking, likely David Fouhey, in front of a whiteboard.

Image + Text

MDNet: A Semantically and Visually Interpretable Medical Image Diagnosis Network

Zizhao Zhang, Yuanpu Xie, Fuyong Xing, Mason McGough, Lin Yang
University of Florida
zizhao@cise.ufl.edu




Figure 1: Overview of our medical image diagnosis network (MDNet) for interpretable diagnosis process.

Image + Text




Figure 2: Overall illustration of MDNet. We use a bladder image with its diagnostic report as an example. The image model generates an image feature to pass to LSTM in the form of a task tuple and a Conv feature embedding (for the attention model) computed by the AAS module (defined in the method). LSTM executes prediction tasks according to the specified image feature type (best viewed in color).

Image + Text




Figure 3: The attention maps of the original method (middle) and our method (right). Our method generates more focal attention on informative (urothelial) regions.






Figure 6: The illustration of class-specific attention. From top to bottom, test images, pathologist annotations, and class attention maps. Like the pathologist annotations, the attention maps are most activated in urothelial regions, largely ignoring stromal or background regions. Best viewed in color.

Strategy



↓ CNN




↓ DNN or RNN

↓ DNN or RNN

Cancer?

MHVZ0ED, TEST	CONFIDENTIAL	Page 11 of 27
<i>VA Medication History</i>		
Source: VA		
Last Updated: 18-Jul 2012 @ 0002		
Sorted By: Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.		
Medication: DOCUSATE NA 100MG CAP		
Instructions: TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status: Active		
Refills Remaining: 10		
Last Filled On: 12-Apr-2012		
Initially Ordered On: 12-Apr-2012		
Quantity	Days Supply	Pharmacy
30	30	PORTLAND PHARMACY
Prescription Number: 11305543A		
Medication: METOPROLOL TARTRATE 50MG TAB		
Instructions: TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status: Active		
Refills Remaining: 3		
Last Filled On: 12-Apr-2012		

Strategy



↓ CNN



↓ DNN or RNN

Segmentation

VA Medication History			
Source:	VA		
Last Updated:	18-Jul 2012 @ 0002		
Sorted By:	Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCUSATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12-Apr-2012		
Initially Ordered On:	12-Apr-2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A
Medication:	METOPROLOL TARTRATE 50MG TAB		
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status:	Active		
Refills Remaining:	3		

CONFIDENTIAL

Page 11 of 27

↓ DNN or RNN

Bimodal Network Architectures for Automatic Generation of Image Annotation from Text

Mehdi Moradi^(✉), Ali Madani, Yaniv Gur, Yufan Guo,
and Tanveer Syeda-Mahmood

IBM Research - Almaden Research Center, San Jose, USA
mmoradi@us.ibm.com

Image
(with annotation marks during training)



Four convolutional layers
[convolution, max-pooling, drop-out]

Two fully connected layers
[512,256]

Medical Report

Stable cardiomedastinal silhouette with mild cardiomegaly.

MeSH Coding


Cardiomegaly/mild

Vector Embedding
(GloVe)

LSTM block

concatenate

Three layer
fully connected
Network
[512, 256, 256]






Fig. 2. Architecture 2: Image and text networks are trained separately. The resulting feature vectors are used to train a network that estimates ROI coordinates.





Fig. 3. Sample results from Architecture 2. Red is the estimated quadrilateral, blue is the one marked by a radiologist.

Strategy



↓ CNN



↓ DNN or RNN

Cancer?


VA Medication History			
Source:	VA		
Last Updated:	18-Jul 2012 @ 0002		
Sorted By:	Last Filled On [Descending]		
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCUSATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12-Apr-2012		
Initially Ordered On:	12-Apr-2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A
Medication:	METOPROLOL TARTRATE 50MG TAB		
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS		
Status:	Active		
Refills Remaining:	3		

CONFIDENTIAL

Page 11 of 27

↓ DNN or RNN

Strategy



CNN, DNN, RNN



VA Medication History			
Source:	VA	Last Updated:	18 Jul 2012 @ 0002
Sorted By:	Last Filled On [Descending]	VA Medication History includes up to two years of medication history unless you select a different date range in your download request.	
<hr/>			
Medication:	DOLUSATE NA 100MG CAP	Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY
Status:	Active	Refills Remaining:	10
Initially Ordered On:	12 Apr 2012	Last Filled On:	12 Apr 2012
Quantity	30	Days Supply	30
		Pharmacy	PORLAND PHARMACY
		Prescription Number	11305543A
<hr/>			
Medication:	METOPROLOL TARTRATE 50MG TAB	Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS
Status:	Active	Refills Remaining:	3
<small>Last Edited On: 12 Apr 2012</small>			

Image to Text

On the Automatic Generation of Medical Imaging Reports

Baoyu Jing^{†*} Pengtao Xie^{†*} Eric P. Xing[†]

[†]Petuum Inc, USA

*School of Computer Science, Carnegie Mellon University, USA

{baoyu.jing, pengtao.xie, eric.xing}@petuum.com



Impression: No acute cardiopulmonary abnormality.

Findings: There are no focal areas of consolidation. No suspicious pulmonary opacities. Heart size within normal limits. No pleural effusions. There is no evidence of pneumothorax. Degenerative changes of the thoracic spine.

MTI Tags: degenerative change

Image to Text




Figure 2: Illustration of the proposed model. MLC denotes a *multi-label classification* network. Semantic features are the word embeddings of the predicted tags. The boldfaced tags “calcified granuloma” and “granuloma” are attended by the co-attention network.

Image to Text



Ground Truth

No active disease. The heart and lungs have in the interval. Both lungs are clear and expanded. Heart and mediastinum normal.



No evidence of active disease. The lungs are clear. There is no focal airspace consolidation. No pleural effusion or pneumothorax. Heart size and mediastinal contour are within normal limits. There are multilevel degenerative changes of the spine.



No acute cardiopulmonary abnormality. Normal heart size mediastinal contours. Eventration of the right hemidiaphragm. No focal airspace consolidation. No pleural effusion or pneumothorax.

Ours-CoAttention

No active disease. The heart and lungs have in the interval. Lungs are clear and expanded. Cardiomedastinal silhouette is within normal limits. No pleural effusion or pneumothorax is seen. No pleural effusion. No cavitary or pneumothorax.

No acute cardiopulmonary findings. Heart size is not enlarged. No focal airspace consolidation suspicious pulmonary opacity large pleural effusion or pneumothorax. No focal areas of consolidation. Degenerative changes of the spine. This is moderate exam of the hydropneumothorax. Lungs are clear. There is no focal airspace consolidation pleural effusion or pneumothorax.

No acute cardiopulmonary abnormality. Stable appearance of the thoracic aorta. The right lateral lower lobe is noted in the right lower right midlung. No large pleural effusion or focal airspace disease. Mild interstitial opacities. Atherosclerotic calcifications bony structures bilaterally. There is no pleural effusion or pneumothorax developed in the right lower lobe.

Ours-no-Attention

The lungs are clear bilaterally. They are grossly normal. No focal lung consolidation. No acute bony abnormality. cm nodule within the right lower lobe on the lateral view. No pneumothorax or pleural effusion. No acute bony abnormality. The heart is not enlarged. The lungs are clear. No acute bony abnormality.

The lungs are clear bilaterally. They are grossly normal. No pleural effusion. The heart is normal in size and contour. The lungs are clear. There are no acute bony findings.

The lungs are clear bilaterally. They are grossly normal. No acute bony abnormality. The lungs are otherwise clear. No acute osseous abnormality. No pleural effusion or pneumothorax. Heart size and pulmonary vascularity appear within normal limits.

Soft Attention

No acute cardiopulmonary abnormality. The lungs are clear bilaterally. Specifically no evidence of focal airspace consolidation pleural effusion or pneumothorax. Cardio mediastinal silhouette is unremarkable. Visualized osseous structures of the thorax are without acute abnormality.

No acute cardiopulmonary abnormality. The lungs are clear bilaterally. There is no pleural effusion or pneumothorax. The heart and mediastinum are normal. There is no focal air space opacity to suggest a pneumonia.

No acute cardiopulmonary abnormality. The lungs are clear bilaterally. There is no focal airspace consolidation. No pleural effusion or pneumothorax. Heart size and pulmonary vascularity appear within normal limits.

Image to Text

TieNet: Text-Image Embedding Network for Common Thorax Disease Classification and Reporting in Chest X-rays

Xiaosong Wang^{*1}, Yifan Peng^{*2}, Le Lu¹, Zhiyong Lu², Ronald M. Summers¹

¹Department of Radiology and Imaging Sciences, Clinical Center,

² National Center for Biotechnology Information, National Library of Medicine,
National Institutes of Health, Bethesda, MD 20892

{xiaosong.wang,yifan.peng,luzh,rms}@nih.gov, lel@nvidia.com

Towards Automated Reporting for Chest X-ray




Figure 1. Overview of the proposed automated chest X-ray reporting framework. A multi-level attention model is introduced.

Image to Text




Image to Text

		Image Sample cases			
	P	Atelectasis Effusion	No finding	Nodule Pneumothorax Mass Consolidation	Mass
Original report		findings : a single ap view of the chest demonstrates increasing bibasilar interstitial opacities with decreased overall aeration . increasing blunting of right costophrenic angle. impression : increasing bibasilar atelectasis with possible development of right pleural effusion .	Normal no evidence of lung infiltrate .	findings : heart and mediastinum unchanged . multiple lung nodules . evidence of recent left chest surgery with left chest tube in place . very small left apical pneumothorax . lungs unchanged , no evidence of acute infiltrates . impression : stable chest .	findings : large left suprähilar and infrähilar masses as well as the well circumscribed nodule the level of the aortic knob . the right infrähilar mass as well . no effusion . impression : metastatic lung disease .
Generated Report		findings : a single ap view of the chest demonstrates unchanged bilateral reticular opacities , consider atelectasis . continued left basilar atelectasis . no evidence of developing infiltrate . the cardiac and mediastinal contours are stable . impression : no evidence of developing infiltrate .	findings : pa and lateral views of the chest demonstrate lungs that are clear without focal mass , infiltrate or effusion . cardiomedastinal silhouette is normal size and contour . pulmonary vascularity is normal in caliber and distribution . impression : no evidence of acute pulmonary pathology	findings : pa and lateral views of the chest demonstrate unchanged bilateral chest tubes . again pulmonary nodules are seen on the right and cardiac silhouette unchanged . the cardiac and mediastinal contours are stable . impression : 1. bilateral masses and left lower lung field consolidation . 2.new bilateral lung masses .	comparison is to previous upright study of no significant interval change is seen in the appearance of the chest . the mediastinal soft tissue and pulmonary vascularity are stable . there are blastic bone lesions in the chest . bones , soft tissues are normal . the lung fields are clear . there are calcified lymph nodes in the left lower lung . impression : sclerotic lesions in the left humeral , consistent with metastasis.

Image to Text

Multimodal Recurrent Model with Attention for Automated Radiology Report Generation

Yuan Xue¹, Tao Xu², L. Rodney Long³, Zhiyun Xue³, Sameer Antani³,
George R. Thoma³, and Xiaolei Huang¹(✉)

¹ College of Information Sciences and Technology,
Penn State University, University Park, PA, USA

sharon.x.huang@gmail.com


² Department of Computer Science and Engineering,
Lehigh University, Bethlehem, PA, USA

³ National Library of Medicine, National Institutes of Health, Bethesda, MD, USA

Abstract. Radiologists routinely examine medical images such as X-Ray, CT, or MRI and write reports summarizing their descriptive findings and conclusive impressions. A computer-aided radiology report generation system can lighten the workload for radiologists considerably and assist them in decision making. Although the rapid development of deep learning technology makes the generation of a single conclusive sentence possible, results produced by existing methods are not sufficiently reliable due to the complexity of medical images. Furthermore, generating detailed paragraph descriptions for medical images remains a challenging problem. To tackle this problem, we propose a novel generative model which generates a complete radiology report automatically. The proposed model incorporates the Convolutional Neural Networks (CNNs) with the Long Short-Term Memory (LSTM) in a recurrent way. It is capable of not only generating high-level conclusive impressions, but also generating detailed descriptive findings sentence by sentence to support the conclusion. Furthermore, our multimodal model combines the encoding of the image and one generated sentence to construct an attention input to guide the generation of the next sentence, and henceforth maintains coherence among generated sentences. Experimental results on the publicly available Indiana U. Chest X-rays from the Open-i image collection show that our proposed recurrent attention model achieves significant improvements over baseline models according to multiple evaluation metrics.

Image to Text

Input Image



Recurrent Attention

Findings: The heart size and mediastinal contours appear within normal limits. No focal airspace consolidation , pleural effusion or pneumothorax. No acute bony abnormalities.

Impression: No acute cardiopulmonary finding.

Findings: The heart size and mediastinal silhouette are within normal limits for contour. The lungs are clear. No focal airspace consolidation. No pleural effusion or pneumothorax. Normal cardiomedastinal silhouette. Heart size is normal.

Impression: Clear lungs. No acute cardiopulmonary abnormality.

Ground Truth

Findings: The heart size and mediastinal silhouette are within normal limits for contour. The lungs are clear. No pneumothorax or pleural effusions. The XXXX are intact.

Impression: No acute cardiopulmonary abnormalities.

Findings: Mediastinal contours are within normal limits. Heart size is within normal limits. No focal consolidation, pneumothorax or pleural effusion. No bony abnormality. Vague density in right mid lung, XXXX related to scapular tip and superimposed ribs. Not visualized on lateral exam.

Impression: Vague density in right XXXX, XXXX related to scapular tip and superimposed ribs. Consider oblique images to exclude true nodule. 2. No acute cardiopulmonary abnormality.

Fig. 1. Examples of original reports vs. reports generated by our recurrent attention model. Note that, Findings is a paragraph containing some descriptive sentences; Impression is a conclusive sentence. XXXXs are wrongly removed keywords due to de-identification.

Image to Text





Fig. 2. The architecture of the proposed multimodal recurrent generation model with attention for radiology reports. Best viewed in color.


Strategy



CNN, DNN, RNN

VA Medication History			
Source: VA			
Last Updated: 18 Jul 2012 @ 0002			
Sorted By: Last Filled On [Descending]			
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.			
Medication:	DOCUSATE NA 100MG CAP	Days Supply:	30
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11305543A
Refills Remaining:	10	Last Filled On:	12 Apr 2012
Initially Ordered On:	12 Apr 2012		
Quantity	30	Days Supply	30
		Pharmacy	PORTLAND PHARMACY
		Prescription Number	11305543A
Medication:	METOPROLOL TARTRATE 50MG TAB	Days Supply:	30
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS	Pharmacy:	PORTLAND PHARMACY
Status:	Active	Prescription Number:	11305543A
Refills Remaining:	3	Last Filled On:	12 Apr 2012
Initially Ordered On:	12 Apr 2012		
Quantity	30	Days Supply	30
		Pharmacy	PORTLAND PHARMACY
		Prescription Number	11305543A

Strategy



CNN, DNN, RNN



CONFIDENTIAL

MHVZZOED , TEST

VA Medication History

Source: VA
Last Updated: 18 Jul 2012 @ 0002
Sorted By: Last Filled On (Descending)
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.

Medication:	DOCUSTATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12 Apr 2012		
Initially Ordered On:	12 Apr 2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A

Medication:	METOPROLOL TARTRATE 50MG TAB
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS
Status:	Active
Refills Remaining:	3
Last Filled On:	12 Apr 2012



CONFIDENTIAL

MHVZZOED , TEST

VA Medication History

Source: VA
Last Updated: 18 Jul 2012 @ 0002
Sorted By: Last Filled On (Descending)
VA Medication History includes up to two years of medication history unless you select a different date range in your download request.


Medication:	DOCUSTATE NA 100MG CAP		
Instructions:	TAKE ONE CAPSULE BY MOUTH EVERY DAY		
Status:	Active		
Refills Remaining:	10		
Last Filled On:	12 Apr 2012		
Initially Ordered On:	12 Apr 2012		
Quantity	Days Supply	Pharmacy	Prescription Number
30	30	PORTLAND PHARMACY	11305543A

Medication:	METOPROLOL TARTRATE 50MG TAB
Instructions:	TAKE ONE TABLET BY MOUTH EVERY 12 HOURS
Status:	Active
Refills Remaining:	3
Last Filled On:	12 Apr 2012

CNN, DNN, RNN




GAN



Generator



True



Discriminator

<https://arxiv.org/abs/1511.06434>

Text to Image

Generative Adversarial Text to Image Synthesis

**Scott Reed, Zeynep Akata, Xinchen Yan, Lajanugen Logeswaran
Bernt Schiele, Honglak Lee**

¹ University of Michigan, Ann Arbor, MI, USA (UMICH.EDU)

this small bird has a pink breast and crown, and black primaries and secondaries.



this magnificent fellow is almost all black with a red crest, and white cheek patch.



the flower has petals that are bright pinkish purple with white stigma



this white and yellow flower have thin white petals and a round yellow stamen



Text to Image




Figure 2. Our text-conditional convolutional GAN architecture. Text encoding $\varphi(t)$ is used by both generator and discriminator. It is projected to a lower-dimensions and depth concatenated with image feature maps for further stages of convolutional processing.

Text to Image

	GT	Ours		GT	Ours		GT	Ours
a group of people on skis stand on the snow.			a man in a wet suit riding a surfboard on a wave.			a pitcher is about to throw the ball to the batter.		
a table with many plates of food and drinks			two plates of food that include beans, guacamole and rice.			a picture of a very clean living room.		
two giraffe standing next to each other in a forest.			a green plant that is growing out of the ground.			a sheep standing in a open grass field.		
a large blue octopus kite flies above the people having fun at the beach.			there is only one horse in the grassy field.			a toilet in a small room with a window and unfinished walls.		

Text to Image

AttnGAN: Fine-Grained Text to Image Generation with Attentional Generative Adversarial Networks

Tao Xu^{*1}, Pengchuan Zhang², Qiuyuan Huang²,
Han Zhang³, Zhe Gan⁴, Xiaolei Huang¹, Xiaodong He⁵

¹Lehigh University ²Microsoft Research ³Rutgers University ⁴Duke University ⁵JD AI Research
{tax313, xih206}@lehigh.edu, {penzhan, qihua, xiaohe}@microsoft.com
han.zhang@cs.rutgers.edu, zhe.gan@duke.edu, xiaodong.he@jd.com

Text to Image





Figure 2. The architecture of the proposed AttnGAN. Each attention model automatically retrieves the conditions (*i.e.*, the most relevant word vectors) for generating different sub-regions of the image; the DAMSM provides the fine-grained image-text matching loss for the generative network.

Text to Image


the bird has a **yellow** crown and a **black** eyering that **is round**



1:bird 4:yellow 0:the 12:round 11:is




1:bird 4:yellow 0:the 8:black 12:round



this bird has a **green** crown **black** primaries and a **white** belly



1:bird 0:this 2:has 11:belly 10:white



6:black 4:green 10:white 0:this 1:bird



Text to Image

<https://t2i.cvalenzuelab.com/>

Integration of Spatial Distribution in Imaging-Genetics

Vaishnavi Subramanian^{1(✉)}, Weizhao Tang¹, Benjamin Chidester², Jian Ma²,
and Minh N. Do¹

¹ Department of Electrical and Computer Engineering, University of Illinois at Urbana-Champaign, Urbana, IL, USA
`vs5@illinois.edu`

² School of Computer Science, Carnegie Mellon University, Pittsburgh, PA, USA

Abstract. To better understand diseases such as cancer, it is crucial for computational inference to quantify the spatial distribution of various cell types within a tumor. To this end, we used Ripley’s K-statistic, which captures the spatial distribution patterns at different scales of — both individual point sets and interactions between multiple point sets. We propose to improve the expressivity of histopathology image features by incorporating this descriptor to capture potential cellular interactions, especially interactions between lymphocytes and epithelial cells. We demonstrate the utility of the Ripley’s K-statistic by analyzing digital slides from 710 TCGA breast invasive carcinoma (BRCA) patients. In particular, we consider its use in the context of imaging-genetics to understand correlations between gene expression and image features using canonical correlation analysis (CCA). Our analysis shows that including these spatial features leads to more significant associations between image features and gene expression.







Fig. 1. (a) Our work-flow comprising CNN-based image feature extraction, spatial feature descriptor, and CCA between the features. (b) Pictorial representation of the K-function evaluated at radius r for the blue process, by counting blue points (self K-function) and red stars (cross K-function) within radius r .




(a) Configuration 1



(b) Configuration 2



(c) Corresponding K-functions



(a) Correlation of spatial-augmented image features with variates (K-function evaluated for $r \leq 500$). (b) Ordered correlations of gene with spatial variate.

Fig. 3. SCCA for gene expression and image features, with spatial features added. (a) The fourth variate has a strong negative correlation with most spatial features. (b) The correlation of gene expression with this variate showed very high correlation with a few genes and then a linear decay in correlation. We chose the transition at the 480th gene, corresponding to a correlation threshold of 0.4, for use in pathway analysis.

TextRay: Mining Clinical Reports to Gain a Broad Understanding of Chest X-Rays

Jonathan Laserson^{1(✉)}, Christine Dan Lantsman², Michal Cohen-Sfady¹, Itamar Tamir³, Eli Goz¹, Chen Brestel¹, Shir Bar⁴, Maya Atar⁵, and Eldad Elnekave¹

¹ Zebra Medical Vision LTD, Shefayim, Israel

jonil@zebra-med.com


² Sheba Medical Center and Tel Aviv University, Ramat Gan, Israel

³ Rabin Medical Center, Petah Tikva, Israel

⁴ Technion, Israel Institute of Technology, Haifa, Israel

⁵ Ben Gurion University, Beersheba, Israel

Abstract. The chest X-ray (CXR) is by far the most commonly performed radiological examination for screening and diagnosis of many cardiac and pulmonary diseases. There is an immense world-wide shortage of physicians capable of providing rapid and accurate interpretation of this study. A radiologist-driven analysis of over two million CXR reports generated an ontology including the 40 most prevalent pathologies on CXR. By manually tagging a relatively small set of sentences, we were able to construct a training set of 959k studies. A deep learning model was trained to predict the findings given the patient frontal and lateral scans. For 12 of the findings we compare the model performance against a team of radiologists and show that in most cases the radiologists agree on average more with the algorithm than with each other.



Training Labels

0	cardiomegaly
0	pulmonary edema
1	vertebral height loss
0	consolidation
0	cardiac pacer
0	rib fracture
	...
0	pleural effusion
0	pneumothorax
1	hilar
0	atelectasis
1	granuloma

Heart size is normal. Mediastinal width is within normal limits. No edema. No focal infiltrate. No pleural effusion or pneumothorax. Right hilar and right lung base calcifications. There is a very mild anterior wedge deformity of a midthoracic vertebrae, possibly T7. Correlate for midthoracic tenderness. No displaced, acute rib fractures are identified.

Fig. 1. TextRay Model Illustration. Frontal (PA) and lateral view images each go through a separate CNN. A fully-connected layer is applied on their concatenated feature vectors and emits the confidence for each finding. Training labels were extracted by analyzing the report sentences. Negative (green) and positive (red) sentences identified. Findings in positive sentences receive a positive training label. Negative or unmentioned findings receive a negative label.

Finding	Pool		Avg. agreement w/ rads			Δ (CI)
	pos	neg	Report	Avg. rad	Textray	Textray vs. rads
Pulmonary edema	128	482	0.613	0.639	0.730	+0.09 (0.07, 0.11)
Elevated diaphragm	202	77	0.731	0.675	0.754	+0.08 (0.05, 0.10)
Abnormal aorta	198	80	0.736	0.693	0.771	+0.08 (0.05, 0.11)
Hyperinflation	95	80	0.678	0.619	0.657	+0.04 (-0.02, 0.10)
Vertebral height loss	126	55	0.781	0.742	0.757	+0.02 (-0.02, 0.06)
Atelectasis	201	78	0.778	0.756	0.767	+0.01 (-0.03, 0.04)
Cardiomegaly	238	372	0.755	0.861	0.866	+0.01 (-0.02, 0.03)
Pleural effusion	207	73	0.905	0.893	0.896	+0.00 (-0.02, 0.03)
Consolidation	194	78	0.690	0.730	0.707	-0.02 (-0.07, 0.02)
Pneumothorax	111	124	0.830	0.855	0.823	-0.03 (-0.08, 0.01)
Rib fracture	183	76	0.683	0.799	0.745	-0.05 (-0.10, -0.01)
Hilar prominence	184	426	0.552	0.797	0.736	-0.06 (-0.09, -0.03)

Table 2 shows that TextRay is on par with human radiologists (within the 95% CI) on 10 out of 12 findings, with the exception of *rib fracture* and *hilar prominence*. On some findings (*elevated diaphragm*, *abnormal aorta*, and *pulmonary edema*), radiologists agree significantly more with our algorithm than

Ordinal Multi-modal Feature Selection for Survival Analysis of Early-Stage Renal Cancer

Wei Shao¹, Jun Cheng², Liang Sun¹, Zhi Han⁴, Qianjin Feng³,
Daoqiang Zhang^{1(✉)}, and Kun Huang^{4(✉)}

College of Computer Science and Technology, Nanjing University of Aeronautics
and Astronautics, Nanjing 211106, China
dqzhang@nuaa.edu.cn

School of Biomedical Engineering, Shenzhen University, Shenzhen 518073, China

³ School of Biomedical Engineering, Southern Medical University, Guangzhou
510515, China

⁴ School of Medicine, Indiana University, Indianapolis, IN 46202, USA
kunhuang@iu.edu

Abstract. Existing studies have demonstrated that combining genomic data and histopathological images can better stratify cancer patients with distinct prognosis than using single biomarker, for different biomarkers may provide complementary information. However, these multi-modal data, most high-dimensional, may contain redundant features that will deteriorate the performance of the prognosis model, and therefore it has become a challenging problem to select the informative features for survival analysis from the redundant and heterogeneous feature groups. Existing feature selection methods assume that the survival information of one patient is independent to another, and thus miss the ordinal relationship among the survival time of different patients. To solve this issue, we make use of the important ordinal survival information among different patients and propose an ordinal sparse canonical correlation analysis (i.e., OSCCA) framework to simultaneously identify important image features and eigengenes for survival analysis. Specifically, we formulate our framework basing on sparse canonical correlation analysis model, which aims at finding the best linear projections so that the highest correlation between the selected image features and eigengenes can be achieved. In addition, we also add constrains to ensure that the ordinal survival information of different patients is preserved after projection. We evaluate the effectiveness of our method on an early-stage renal cell carcinoma dataset. Experimental results demonstrate that the selected features correlated strongly with survival, by which we can achieve better patient stratification than the comparing methods.




Fig. 1. The flowchart of the proposed method.




Fig. 2. Comparisons of the survival curves by applying different feature methods.

DeepEM: Deep 3D ConvNets with EM for Weakly Supervised Pulmonary Nodule Detection

Wentao Zhu¹(✉), Yeeleng S. Vang¹, Yufang Huang², and Xiaohui Xie¹

¹ University of California, Irvine, USA

{wentaoz1,ysvang,xhx}@uci.edu

² Lenovo Research, Beijing, China
huangyf8@lenovo.com

Abstract. Recently deep learning has been witnessing widespread adoption in various medical image applications. However, training complex deep neural nets requires large-scale datasets labeled with ground truth, which are often unavailable in many medical image domains. For instance, to train a deep neural net to detect pulmonary nodules in lung computed tomography (CT) images, current practice is to manually label nodule locations and sizes in many CT images to construct a sufficiently large training dataset, which is costly and difficult to scale. On the other hand, electronic medical records (EMR) contain plenty of partial information on the content of each medical image. In this work, we explore how to tap this vast, but currently unexplored data source to improve pulmonary nodule detection. We propose DeepEM, a novel deep 3D ConvNet framework augmented with expectation-maximization (EM), to mine weakly supervised labels in EMRs for pulmonary nodule detection. Experimental results show that DeepEM can lead to 1.5% and 3.9% average improvement in free-response receiver operating characteristic (FROC) scores on LUNA16 and Tianchi datasets, respectively, demonstrating the utility of incomplete information in EMRs for improving deep learning algorithms (<https://github.com/uci-cbcl/DeepEM-for-Weakly-Supervised-Detection.git>).




Fig. 1. Illustration of DeepEM framework. Faster R-CNN is employed for nodule proposal generation. Half-Gaussian model and logistic regression are employed for central slice and lobe location respectively. In the E-step, we utilize all the observations, CT slices, and weak label to infer the latent variable, nodule proposals, by maximum a posteriori (MAP) or sampling. In the M-step, we employ the estimated proposals to update parameters in the Faster R-CNN and logistic regression.



















Ground truth center slice		slice 83 center (48,240) diam. 9		slice 114 center (40,88) diam. 7		slice 129 center (129,99) diam. 14		slice 114 center (137,75) diam. 29		slice 47 center (122,202) diam. 12
Faster R-CNN center slice		slice 84 center (48,242) diam. 10		slice 113 center (41,89) diam. 9		slice 126 center (127,98) diam. 19		slice 111 center (138,74) diam. 36		slice 47 center (122,202) diam. 17
DeepEM (MAP) center slice		slice 84 center (48,241) diam. 8		slice 113 center (41,88) diam. 5		slice 128 center (129,98) diam. 16		slice 113 center (137,73) diam. 28		slice 48 center (122,202) diam. 15
DeepEM (sampling) center slice		slice 84 center (48,241) diam. 8		slice 113 center (41,88) diam. 5		slice 128 center (129,98) diam. 16		slice 113 center (137,73) diam. 29		slice 46 center (121,201) diam. 13

Fig. 3. Detection visual comparison among Faster R-CNN [14], DeepEM with MAP and DeepEM with Sampling on nodules randomly sampled from Tianchi. DeepEM provides more accurate detection (central slice, center and diameter) than Faster R-CNN.

Towards Automatic Report Generation in Spine Radiology Using Weakly Supervised Framework

Zhongyi Han^{1,2}, Benzheng Wei^{1,2(✉)}, Stephanie Leung^{3,4}, Jonathan Chung^{3,4}, and Shuo Li^{3,4(✉)}

¹ College of Science and Technology, Shandong University of Traditional Chinese Medicine, Jinan, Shandong, China
wbz99@sina.com

² Computational Medicine Lab (CML), Shandong University of Traditional Chinese Medicine, Jinan, Shandong, China

³ Department of Medical Imaging, Western University, London, ON, Canada
slishuo@gmail.com

⁴ Digital Imaging Group (DIG), London, ON, Canada

Abstract. The objective of this work is to automatically generate unified reports of lumbar spinal MRIs in the field of radiology, i.e., given an MRI of a lumbar spine, directly generate a radiologist-level report to support clinical decision making. We show that this can be achieved via a weakly supervised framework that combines deep learning and symbolic program synthesis theory to overcome four inevitable tasks: semantic segmentation, radiological classification, positional labeling, and structural captioning. The weakly supervised framework using object level annotations without requiring radiologist-level report annotations to generate unified reports. Each generated report covers almost type lumbar structures comprised of six intervertebral discs, six neural foramina, and five lumbar vertebrae. The contents of each report contain the exact locations and pathological correlations of these lumbar structures as well as their normalities in terms of three type relevant spinal diseases: intervertebral disc degeneration, neural foraminal stenosis, and lumbar vertebrae deformities. This framework is applied to a large corpus of T1/T2-weighted sagittal MRIs of 253 subjects acquired from multiple vendors. Extensive experiments demonstrate that the framework is able to generate unified radiological reports, which reveals its effectiveness and potential as a clinical tool to relieve spinal radiologists from laborious workloads to a certain extent, such that contributes to relevant time savings and expedites the initiation of many specific therapies.






Fig. 1. The workflow of the proposed weakly supervised framework.



(a) Radiological Report

At T12-L1, the intervertebral disc has obvious degenerative changes. The neural foramen does not have stenosis.

At L1-L2, the above vertebra has deformative changes. The intervertebral disc does not have obvious degenerative changes. The neural foramen does not have obvious stenosis.

At L2-L3, the neural foramen has obvious stenosis. The intervertebral disc does not have obvious degenerative changes. The above vertebra does not have deformative changes.

At L3-L4, disc degenerative changes are associated with neural foraminal stenosis.

At L4-L5, disc degenerative changes are associated with neural foraminal stenosis.

At L5-S1, the intervertebral disc has obvious degenerative changes. The above vertebra also has deformative changes. They lead to the neural foraminal stenosis.

(b) Radiological Report

At T12-L1, the intervertebral disc has obvious degenerative changes. The neural foramen does not have stenosis.

At L1-L2, the above vertebra does not have deformative changes. The intervertebral disc does not have degenerative changes. The neural foramen also does not have stenosis.

At L2-L3, the above vertebra does not have deformative changes. The intervertebral disc does not have degenerative changes. The neural foramen also does not have stenosis.

At L3-L4, the neural foramen has obvious stenosis. The intervertebral disc does not have obvious degenerative changes. The above vertebra does not have deformative changes.

At L4-L5, disc degenerative changes are associated with neural foraminal stenosis.

At L5-S1, the intervertebral disc has obvious degenerative changes. The above vertebra also has deformative changes. They lead to neural foraminal stenosis to a certain extent.

Fig. 2. The generated radiological reports. The text in purple color represents that our framework is helpful for building comprehensive pathological analysis.

Accurate Weakly-Supervised Deep Lesion Segmentation Using Large-Scale Clinical Annotations: Slice-Propagated 3D Mask Generation from 2D RECIST

Jinzheng Cai^{1,2}(✉), Youbao Tang¹, Le Lu¹, Adam P. Harrison¹, Ke Yan¹,
Jing Xiao³, Lin Yang², and Ronald M. Summers¹

¹ National Institutes of Health Clinical Center, Bethesda, MD 20892, USA

{youbao.tang, le.lu, adam.harrison, ke.yan, rms}@nih.gov

² University of Florida, Gainesville, FL 32611, USA
jimmycai@ufl.edu, lin.yang@bme.ufl.edu

³ Ping An Insurance (Group) Company of China, Ltd.,
Shenzhen 510852, People's Republic of China
xiaojing661@pingan.com.cn

Abstract. Volumetric lesion segmentation from computed tomography (CT) images is a powerful means to precisely assess multiple time-point lesion/tumor changes. However, because manual 3D segmentation is prohibitively time consuming, current practices rely on an imprecise surrogate called response evaluation criteria in solid tumors (RECIST). Despite their coarseness, RECIST markers are commonly found in current hospital picture and archiving systems (PACS), meaning they can provide a potentially powerful, yet extraordinarily challenging, source of weak supervision for full 3D segmentation. Toward this end, we introduce a convolutional neural network (CNN) based weakly supervised slice-propagated segmentation (WSSS) method to (1) generate the initial lesion segmentation on the axial RECIST-slice; (2) learn the data distribution on RECIST-slices; (3) extrapolate to segment the whole lesion slice by slice to finally obtain a volumetric segmentation. To validate the proposed method, we first test its performance on a fully annotated lymph node dataset, where WSSS performs comparably to its fully supervised counterparts. We then test on a comprehensive lesion dataset with 32,735 RECIST marks, where we report a mean Dice score of 92% on RECIST-marked slices and 76% on the entire 3D volumes.




Fig. 1. Overview of the proposed method. **Right:** we use CNN outputs to gradually generate extra training data for lesion segmentation. Arrows colored in red, orange, and blue indicate slice-propagated training at its 1st, 2nd, and 3rd steps, respectively. **Left:** regions colored with red, orange, green, and blue inside the initial segmentation mask Y present FG, PFG, PBG, and BG, respectively.